



INOXI CLEAN

RESTORE
renews corroded surfaces

INOXI CLEAN RESTORE

Good as new!

Safeguard your investment with regular maintenance.

Remove flash rust and other impurities from stainless steel, aluminium & copper in a trice with InoxiClean Restore. Apply InoxiClean Restore; let it work in for 5 minutes, rub and rinse, and everything looks as good as new.

InoxiClean works like a deep cleaner which enables you to quickly remove dirt from deep within the pores of the surface, without affecting the oxide protective layer on the stainless steel surface. Acid rain, soot particles, fine dust, and other unfavorable environmental factors affect the surface of stainless steel and make it lustreless. InoxiClean Restore removes corrosion, dirt and oil, and brings back the material's original lustre in the wink of an eye.

Moreover, InoxiClean Restore applies a water and dirt repellent protective layer, keeping the surface protected from new dirt and corrosion for longer periods of time.

InoxiClean Restore is 100% biodegradable, and the product is therefore not harmful to man or the environment. InoxiClean Restore is suitable for cleaning machines and tools for use in the food industry.

Advantages of InoxiClean Restore:

- Removes rust and impurities in a fast and safe manner
- Brings back the original lustre and appearance of your stainless steel
- 100% biodegradable

Applications:

- Industrial kitchens
- Swimming pools
- Facade renovation
- Street furniture
- Terrace furniture
- Food containers
- Industrial bakeries
- Maintenance of machinery
- Rims and exhaust systems in vehicles



HOW TO USE?

①

Spray the Restore on the surface to be treated.
Leave it to work in approx. 5 minutes.



②

Rub down, applying some pressure,
using a white tex pad.



③

Remove the product with a wet cloth and
subsequently dry using a clean cloth.





RESTORE SET

With the starter kit you have everything at hand to clean all surfaces comfortable and in a jiffy.

Restore

Description	Code
InoxiClean Restore 500 ml	ICRE/500
InoxiClean Restore 5 l	ICRE/5000
InoxiClean Restore kit 500 ml	ICRE/SET

Content ICRE/SET

Description	Quantity	Code
InoxiClean Restore 500 ml	1	ICRE/500
Tex 00 sheet 140 x 100	8	T/00/140100
Tex handle	1	HG12080
microfiber cloth extra fine	1	ICMV/XF/400400
microfiber cloth fine	1	ICMV/F/400400
Gloves (per pair)	3	-



MORE INFORMATION?

Surf to www.cibo.be for more information, product videos, the full range of InoxiClean products and a point of sale in your neighbourhood.



Deugenietaat 5
3150 Tildonk Belgium
T +32 (0)16 61 85 85
F +32 (0)16 61 84 84
info@cibo.be